

# Red Blood cell count result





Entity: Cluster

Concept description:				Identification:	
Test to verify the nr of red blood cells in the body, results as a single value.				<i>Id:</i> openEHR-EHR-CLUSTER.laboratory_test_panel-red_blood_cell_count.v0 <i>Reference model:</i> openEHR_EHR	
Purpose	Use	Misuse	Copyright	References	Contact
To record Red Blood cell count test results as a single value. LOINC 26453-1 Erythrocytes [#/? volume] in Blood Component Property Time System Scale Method Erythrocytes NCnc Pt Bld Qn Part: Erythrocytes Erythrocytes or red blood cells (RBCs) are the cells in the circulation that carry oxygen to and remove carbon dioxide from the tissues throughout the body. They are	To record Red Blood cell count test results as a single value. Normally used in conjunction with a parent Laboratory test result (Observation) archetype.	Should not be used to record Anatomical pathology macroscopic/microscopic findings.	© openEHR Foundation	Based on NEHTA 'Pathology Test' archetype. Available from: <a href="http://dcm.nehta.org.au/ckm/OKM.html#showarchetype_1013.1.839_8">http://dcm.nehta.org.au/ckm/OKM.html#showarchetype_1013.1.839_8</a> Pathology (Data Specifications) Version 1.0 [Internet]. Sydney, Australia: National E-Health Transition Authority; 2007 May 29 [cited 2011 Jul 11]; Available at <a href="http://www.nehta.gov.au/component/docman/doc_download/962-pathology-v10">http://www.nehta.gov.au/component/docman/doc_download/962-pathology-v10</a> . Laboratory Technical Framework, Volume 3: Content, Revision 3.0 [Internet]. USA: IHE International; 2011 May 19; [cited 2011 Jul 11]. Available from: <a href="http://www.ihe.net/Technical_Framework/index.cfm#laboratory_H17">http://www.ihe.net/Technical_Framework/index.cfm#laboratory_H17</a> FHIR Observation resource: HL7 FHIR; Available from <a href="http://www.hl7.org/implement/standards/fhir/observation.html">http://www.hl7.org/implement/standards/fhir/observation.html</a>	

produced in the bone marrow in response to erythropoietin where they transition through six stages over a seven day period. When they are released into the circulation, their nucleus has been extruded and they measure 6-8 microns in diameter. The lifespan of RBCs is about 120 days. When the RBC ages the cell membrane becomes less pliable and tears as they cell travels through the micro vessels of the body. The damaged RBCs are removed from the circulation by the spleen. Variations in number, shape and size of RBCs are indicative of many diseases and disorders. There are many factors that may lead to decreased numbers of RBCs including

hemorrhage,  
hemolysis, iron  
or vitamin  
deficiency,  
marrow failure  
and more. Larger  
than normal  
RBCs may be  
indicative of liver  
disease while  
smaller than  
normal RBCs are  
seen in  
thalassemias and  
other anemias.  
(Mosby's manual  
of diagnostic and  
laboratory tests,  
Kathleen Deska  
Pagana; Timothy  
James Pagana,  
Elsevier St.  
Louis, Mo  
©2010) Source:  
Regenstrief  
LOINC NB: This  
is not cloned in  
templates from  
laboratory-tests  
but specialized.  
The reason for  
this is that  
maintainability  
becomes hard  
when there are  
changes in the  
model, but the  
corresponding  
constraint can  
occur in more  
templates (which  
will happen

because, ET decided to have archetypes which can occur in several templates.

Concept	Description	Constraints	Values
 <b>Laboratory result</b>	Specific detailed result, including both the value of the result item, and additional information that may be useful for clinical interpretation.	<b>Cluster</b> 0..*	
<b>Q</b> <b>Result value</b>	Actual value of the result.	<b>Quantity</b> 0..1	Property = Concentration Units = x 10 <sup>12</sup> /l; >=3.8; <=6.5;
<b>T</b> <b>Comment</b>	Comment about the Result.	<b>Text</b> 0..1	Text;
<b>T</b> <b>Reference range guidance</b>	Additional advice on the applicability of the reference range.	<b>Text</b> 0..1	Text;
<b>T</b> <b>Result status</b>	The status of the result value.	<b>Text</b> 0..1	Internal; 'Registered', 'Interim', 'Final', 'Amended', 'Cancelled/Aborted', 'Not requested'
 <b>Result status timestamp</b>	The date and/or time that the entire result was issued for the recorded 'Result status'.	<b>DateTime</b> 0..1	Allow all
	Slot Result detail [Cluster]	Include : Cluster	Exclude : Cluster
	Slot Other detail [Cluster]	Include : Cluster	Exclude : Cluster