





Sodium (NA+) test result

Entity: Cluster

Concept description:				Identification:	
Sodium (NA+) test results as a single value.				<i>Id:</i> openEHR-EHR-CLUSTER.laboratory_test_panel-sodium.v0 <i>Reference model:</i> openEHR_EHR	
Purpose	Use	Misuse	Copyright	References	Contact
To record Sodium (NA+) test results as a single value. LOINC 2951-2 Sodium [Moles/? volume] in Serum or Plasma Component Sodium Property Time System Scale Method SCnc Pt Ser/Plas Qn Sodium is an essential nutrient that regulates blood volume, blood pressure, osmotic equilibrium and electrolyte balance. Sodium chloride is the principal source of sodium in the diet, and is used for seasoning and	To record Sodium (NA+) test results as a single value. Normally used in conjunction with a parent Laboratory test result (Observation) archetype.	Should not be used to record Anatomical pathology macroscopic/microscopic findings.	© openEHR Foundation	Based on NEHTA 'Pathology Test' archetype. Available from: http://dcm.nehta.org.au/ckm/OKM.html#showarchetype_1013.1.839_8 Pathology (Data Specifications) Version 1.0 [Internet]. Sydney, Australia: National E-Health Transition Authority; 2007 May 29 [cited 2011 Jul 11]; Available at http://www.nehta.gov.au/component/docman/doc_download/962-pathology-v10 . Laboratory Technical Framework, Volume 3: Content, Revision 3.0 [Internet]. USA: IHE International; 2011 May 19; [cited 2011 Jul 11]. Available from: http://www.ihe.net/Technical_Framework/index.cfm#laboratory_H17 FHIR Observation resource: HL7 FHIR; Available from http://www.hl7.org/implement/standards/fhir/observation.html	

as a preservative.
Increased levels
of sodium intake
can cause
hypertension and
reportedly leads
to 7.6 million
premature deaths
worldwide.
Sodium is also
important in
neuron function
and
osmoregulation
between cells and
the extracellular
fluid. NB the
mass/volume unit
is not in LOINC
SNOMED
104934005 |
Sodium
measurement,
serum
(procedure) | NB:
This is not cloned
in templates from
laboratory-tests
but specialized.
The reason for
this is that
maintainability
becomes hard
when there are
changes in the
model, but the
corresponding
constraint can
occur in more
templates (which
will happen
because, ET

decided to have archetypes which can occur in several templates.

Concept	Description	Constraints	Values
 Laboratory result	Specific detailed result, including both the value of the result item, and additional information that may be useful for clinical interpretation.	Cluster 0..*	
Q Result value	Sodium (NA+) test: Actual value of the result.	Quantity 0..1	Property = Concentration Units = mmol/l; >=130; <=155; Units = mg/dl; Units = umol/l;
T Comment	Comment about the Result.	Text 0..1	Text;
T Reference range guidance	Additional advice on the applicability of the reference range.	Text 0..1	Text;
T Result status	The status of the result value.	Text 0..1	Internal; 'Registered', 'Interim', 'Final', 'Amended', 'Cancelled/Aborted', 'Not requested'
 Result status timestamp	The date and/or time that the entire result was issued for the recorded 'Result status'.	DateTime 0..1	Allow all
	Slot Result detail [Cluster]	Include : Cluster	Exclude : Cluster
	Slot Other detail [Cluster]	Include : Cluster	Exclude : Cluster