





# Serum Albumin level

Entity: Cluster

Concept description:				Identification:	
Serum Albumin level test result as a single value.				<i>Id:</i> openEHR-EHR-CLUSTER.laboratory_test_panel-albumin.v0 <i>Reference model:</i> openEHR_EHR	
Purpose	Use	Misuse	Copyright	References	Contact
<p>To record Serum Albumin level test results as a single value. LOINC 1756-6</p> <p>Albumin in CSF/Albumin in Serum or Plasma Component Property Time System Scale Method</p> <p>Albumin.CSF/Albumin.SerPl RelRto Pt Ser/Plas+CSF Qn</p> <p>The albumin quotient (AQ; CSF albumin/serum or plasma albumin x 100) is an important measurement to quantify the blood-brain barrier function. Normal function of the AQ barrier depends on the age of the patient; in newborns, a higher permeability may be observed. Source: Regenstrief LOINC SNOMED 26758005   Albumin measurement (procedure)   NB: This is not cloned in templates from laboratory-tests but specialized. The reason for this is that maintainability becomes hard when there are changes in the model, but the</p>	<p>To record Serum Albumin level test results as a single value. Normally used in conjunction with a parent Laboratory test result (Observation) archetype.</p>	<p>Should not be used to record Anatomical pathology macroscopic/microscopic findings.</p>	<p>© openEHR Foundation</p>	<p>Based on NEHTA 'Pathology Test' archetype. Available from: <a href="http://dcm.nehta.org.au/ckm/OKM.html#showarchetype_1013.1.839_8">http://dcm.nehta.org.au/ckm/OKM.html#showarchetype_1013.1.839_8</a></p> <p>Pathology (Data Specifications) Version 1.0 [Internet]. Sydney, Australia: National E-Health Transition Authority; 2007 May 29 [cited 2011 Jul 11]; Available at <a href="http://www.nehta.gov.au/component/docman/doc_download/962-pathology-v10">http://www.nehta.gov.au/component/docman/doc_download/962-pathology-v10</a>. Laboratory Technical Framework, Volume 3: Content, Revision 3.0 [Internet]. USA: IHE International; 2011 May 19; [cited 2011 Jul 11]. Available from: <a href="http://www.ihe.net/Technical_Framework/index.cfm#laboratory">http://www.ihe.net/Technical_Framework/index.cfm#laboratory</a></p> <p>HL7 FHIR Observation resource: HL7 FHIR; Available from <a href="http://www.hl7.org/implement/standards/fhir/observation.html">http://www.hl7.org/implement/standards/fhir/observation.html</a></p>	<p>Bert Verhees, ROSA Software</p>

corresponding constraint can occur in more templates (which will happen because, ET decided to have archetypes which can occur in several templates.

Concept	Description	Constraints	Values
 <b>Laboratory result</b>	Specific detailed result, including both the value of the result item, and additional information that may be useful for clinical interpretation.	<b>Cluster</b> 0..*	
<b>Q Result value</b>	Actual value of the result.	<b>Quantity</b> 0..1	Property = Concentration Units = mg/dl; >=2500; <=4000; Units = g/l; >=25; <=40;
<b>T Comment</b>	Comment about the Result.	<b>Text</b> 0..1	Text;
<b>T Reference range guidance</b>	Additional advice on the applicability of the reference range.	<b>Text</b> 0..1	Text;
<b>T Result status</b>	The status of the result value.	<b>Text</b> 0..1	Internal; 'Registered', 'Interim', 'Final', 'Amended', 'Cancelled/Aborted', 'Not requested'
 <b>Result status timestamp</b>	The date and/or time that the entire result was issued for the recorded 'Result status'.	<b>DateTime</b> 0..1	Allow all
 <b>Slot Result detail [Cluster]</b>	Slot Result detail [Cluster]	Include : Cluster	Exclude : Cluster
 <b>Slot Other detail [Cluster]</b>	Slot Other detail [Cluster]	Include : Cluster	Exclude : Cluster